



 **GIONO®**  
[www.giono.com.my](http://www.giono.com.my)

**NITRILE GLOVES**

*Your reliable healthcare products*





## GIONO Nitrile Gloves

Giono Nitrile tends to be stronger and more puncture resistant than latex gloves. This essentially means that they are well recommended for patients whose the allergy statuses have not been known. They are tested with similar standards as latex gloves.

If you want to find out about the benefits of using nitrile, do check out our Disposable Gloves Material Guide for a comprehensive breakdown of its standards.

If you have decided to switch from latex to nitrile, but is still unsure of where to start, then you are in the right place. Our guidelines consider beyond the strength, grip, and the durability of gloves. Hopefully, this will give you a better idea of the different types of nitrile gloves that are available for use in a medical environment.



# Product Specifications

## Nitrile Examination Disposable Gloves

Product Description	
Type	Nitrile Examination Glove, Powder Free, Online Single Chlorinated, Non Sterile
Material	100% Synthetic Nitrile Latex
Usage	To conduct medical examination, diagnostic and therapeutic procedures to protect patient and user from cross-contamination or infection
Design and feature	Ambidextrous, beaded cuff, smooth and finger textured
Powder	No powder lubricant added
Color	Blue
Thickness	3.5 to 4.0 mils
Weight of M size	3.5gm (+/- 0.2gm)
Packing	100 pieces per box by weight, 10 inner boxes per carton or 20 inner box per carton
Shelf Life	5 years from the date of manufacture with the above storage condition
Performance standard	EN 455 – Part 1,2 and 3
Certificates	CE, COO and LGM

Description	Size	Standard (mm)	Length (mm)
Palm Width	XS	76 +/- 4	All Size 240 +/-
	S	84 +/- 4	
	M	94 +/- 4	
	L	105 +/- 4	
	XL	113 +/- 4	
Thickness, Single Wall	All Size	Finger: 0.09 +/- 0.02 Typical Value: 0.09 to 0.11	
		Palm: 0.07 +/- 0.02 Typical Value: 0.06 to 0.07	
Inspection	Related Defects	Inspection Level	AQL
Water Tight Test	Holes	G1	1.5
Visual Defects	Major Visual	G1	2.5
Visual Defects	Minor Visual		4.0
Critical Incl	Gloves Counting	S2	4.0

